

PolySyn™/聚乙醇酸

钻石针

钻石针 			规格							
			材质	线长	10-0	9-0	8-0	7-0	6-0	5-0
DE7	8.13mm 3/8 弧		紫色 编织	12"/30cm					G551N△ J551G L1551K	
G-6 HE-6										

铲针

铲针 			规格							
			材质	线长	9-0	8-0	7-0	6-0	5-0	
VSL7 	1/4 弧 针长 6.60mm 曲率 94	半径 3.94mm 弦长 5.89mm 针径 .20mm	紫色 编织	18"/45cm				G544N△ J544G L1750K		
	TG100-8 SE-100-8									
VCL8 	1/4 弧 针长 8.51mm 曲率 112	半径 4.36mm 弦长 7.24mm 针径 .33mm	紫色 编织	8"/20cm				G590N△ J590G L1740K		
	1/4 弧 针长 8.51mm 曲率 109	半径 4.32mm 弦长 7.28mm 针径 .34mm		18"/45cm						G570N J570G ●
DSL7 	3/8 弧 针长 6.6mm 曲率 131	半径 2.79mm 弦长 5.26mm 针径 .2mm	紫色 编织	18"/45cm			GA-546N J546G ●			
				8"/20cm						GA-547N J547G ●
				12"/30cm						GA-548N J548G ●
							GA-566N J566G ●			
DSL6 	3/8 弧 针长 6.15mm 曲率 135	半径 2.54mm 弦长 4.83mm 针径 .15mm	紫色 编织	12"/30cm		GA-974N J974G ●				
DVL04 	3/8 弧 针长 4mm 曲率 145	半径 1.68mm 弦长 3.25mm 针径 .10mm	紫色 单股	5"/13cm	GA-8760N D8760 ●					
VSP8 	1/4 弧 针长 8.51mm 曲率 109	半径 4.32mm 弦长 7.28mm 针径 .34mm	紫色 编织						G553N ● ●	

*蓝色为Sharpoint眼科缝线型号
**黑色及绿色为其他参照型号

M.E.T.™ 针

					规格		
			材质	线长	9-0	8-0	7-0
DTRM4	3/8 弧 针长 4mm	70 Microns 曲率 135	紫色 单股	5"/13cm	GA-402N V402G ●		
DRTM5	3/8 弧 针长 5mm	100 Microns 曲率 135	紫色 编织	5"/13cm		GA-405N V405G ●	
DTRM6	3/8 弧 针长 6mm	140 Microns 曲率 135	紫色 编织	5"/13cm		GA-401N V401G ●	

UltraGlide™超滑铲针

					规格		
			材质	线长	10-0	9-0	8-0
HGL5	1/2 弧 针长 5.5mm 曲率 160	半径 2.03mm 弦长 4.06mm 针径 .15mm	紫色 单股	12"/30cm	GA-448N V448G ●		

*蓝色为Sharp点眼科缝线型号
**黑色及绿色为其他参照型号